



Product Datasheet

Rab 2B antibody (orb770663)


Description	Rabbit polyclonal antibody to Rab 2B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RAB2B. AA range:103-152
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24
Expiration Date	12 months from date of receipt.



Western blot analysis of RAB2B Antibody....



Western blot analysis of the lysates from K562, HeLa, and 293 cells.



Immunohistochemical analysis of paraffin-embedded tissue sections.