



Product Datasheet

TBC1D4 (phospho-Thr642) antibody (orb770994)


Description	Rabbit polyclonal antibody to TBC1D4 (phospho-Thr642)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human AS160 around the phosphorylation site of Thr642. AA range:611-660
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150
Expiration Date	12 months from date of receipt.



Western blot analysis of L929 using p-TB...



Immunofluorescence analysis of HeLa cell...



Immunohistochemistry analysis of paraffi...