

Product Datasheet

OSR1 (phospho-Thr185) antibody (orb771015)

Description Rabbit polyclonal antibody to OSR1 (phospho-Thr185)

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human OSR1 around the phosphorylation site of Thr185. AA range:151-200

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 65

Expiration Date 12 months from date of receipt.



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from He...