

Product Datasheet

Histone H3 (Di Methyl Lys14) antibody (orb771033)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to Histone H3 (Di Methyl Lys14) |
| Species/Host | Rabbit |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Synthetic Peptide of Histone H3 (Di Methyl Lys14) |
| Preservatives | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 15-17 |
| Expiration Date | 12 months from date of receipt. |