

Product Datasheet

Beclin 1 \hat{A} antibody (orb771144)

Description	Rabbit polyclonal antibody to Beclin 1 \hat{A}
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein of Beclin-1
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60
Expiration Date	12 months from date of receipt.