



Product Datasheet

RPA14 antibody (orb771487)

Description	Rabbit polyclonal antibody to RPA14
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide from human protein at AA range: 30-80
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14
Expiration Date	12 months from date of receipt.



Western blot analysis of lysate, antibody...



Immunohistochemical analysis of paraffin...