



Product Datasheet

SYK (Phospho-Tyr348) antibody (orb771750)

Description	Rabbit polyclonal antibody to SYK (Phospho-Tyr348)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized phospho derived from human SYK (Phospho-Tyr348)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	82
Expiration Date	12 months from date of receipt.



Western blot analysis of CACO2 lysate, a...



Immunohistochemical analysis of paraffin...