

Product Datasheet

LATS1/2 (Phospho-Thr1079/1041) antibody (orb771756)

Description Rabbit polyclonal antibody to LATS1/2 (Phospho-Thr1079/1041)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen Synthesized phospho derived from human LATS1/2 (Phospho-Thr1079/1041)

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Isotype IgG

Clonality Polyclonal

MW 140

Expiration Date 12 months from date of receipt.



Western blot analysis of various lysate,...