

Product Datasheet

PTPRB antibody (orb773014)

Description	Rabbit polyclonal antibody to PTPRB
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from human protein . at AA range: 280-360
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	219
Expiration Date	12 months from date of receipt.



Western Blot analysis of HEK293 lysis, u...