

## Product Datasheet

### Mouse ACE ELISA Kit (orb774866)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Mouse
<b>Range</b>	0.32-20 ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<b>P09470</b>
<b>Sensitivity</b>	0.142 ng/mL
<b>Expiration Date</b>	Please enquire.

