

Product Datasheet

Human CRP ELISA Kit (orb774905)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	62.5-4000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, urine, cerebrospinal fluid, cell culture supernates and other biological fluids
Uniprot ID	P02741
Sensitivity	26.1 pg/mL
Expiration Date	Please enquire.

