

Product Datasheet

Mouse C5a ELISA Kit (orb774917)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Mouse |
| Range | 1.57-100 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | P06684 |
| Sensitivity | 0.57 ng/mL |
| Expiration Date | Please enquire. |

