

Product Datasheet

Human LTF ELISA Kit (orb774931)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 0.32-20 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, breast milk, cell culture supernates and other biological fluids |
| Uniprot ID | P02788 |
| Sensitivity | 0.125 ng/mL |
| Expiration Date | Please enquire. |

