

Product Datasheet

Mouse Bid ELISA Kit (orb774932)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P70444
Sensitivity	0.127 ng/mL
Expiration Date	Please enquire.

