

## Product Datasheet

### Canine VEGFA ELISA Kit (orb774942)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Canine
<b>Range</b>	15.63-1000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<b>Q9MYV3</b>
<b>Sensitivity</b>	5.6 pg/mL
<b>Expiration Date</b>	Please enquire.

