

Product Datasheet

Human PDGFAB ELISA Kit (orb774954)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 78.13-5000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Sensitivity | 27 pg/mL |
| Expiration Date | Please enquire. |

