

## Product Datasheet

### Human VEGFR1 ELISA Kit (orb774991)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 78.13-5000 pg/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** serum, plasma, tissue homogenates, saliva, cell lysates, follicular fluid and other biological fluids

**Uniprot ID** **P17948**

**Sensitivity** 29 pg/mL

**Expiration Date** Please enquire.

