

Product Datasheet

Porcine IFN γ ELISA Kit (orb775085)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Porcine
Range	7.82-500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P17803
Sensitivity	2.8 pg/mL
Expiration Date	Please enquire.

