

Product Datasheet

Guinea pig IL1b ELISA Kit (orb775146)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Guinea pig
Range	7.82-500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9WVG1
Sensitivity	2.8 pg/mL
Expiration Date	Please enquire.

