

Product Datasheet

Guinea pig IL8 ELISA Kit (orb775147)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Guinea pig
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P49113
Sensitivity	6.1 pg/mL
Expiration Date	Please enquire.

