

Product Datasheet

Rat ALB ELISA Kit (orb775153)

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| Description | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| Reactivity | Rat |
| Range | 3.13-200 µg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma and other biological fluids |
| Uniprot ID | P02770 |
| Sensitivity | 0.95 µg/mL |
| Expiration Date | Please enquire. |

