

## Product Datasheet

### Gallus Ntn1 ELISA Kit (orb775225)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Gallus
<b>Range</b>	15.63-1000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q90922</a>
<b>Sensitivity</b>	5.4 pg/mL
<b>Expiration Date</b>	Please enquire.

