

Product Datasheet

Mouse FGF21 ELISA Kit (orb775249)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	3.13-200 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9JJN1
Sensitivity	5.2 pg/mL
Expiration Date	Please enquire.

