

Product Datasheet

Porcine INS ELISA Kit (orb775250)

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Porcine |
| Range | 7.82-500 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | P01315 |
| Sensitivity | 2.59 pg/mL |
| Expiration Date | Please enquire. |

