

Product Datasheet

Human ALCAM ELISA Kit (orb775261)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 0.16-10 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | Q13740 |
| Sensitivity | 0.068 ng/mL |
| Expiration Date | Please enquire. |

