

Product Datasheet

Human HSPb2 ELISA Kit (orb775354)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.63-40 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q16082
Sensitivity	0.265 ng/mL
Expiration Date	Please enquire.

