

Product Datasheet

Porcine HSPA1A ELISA Kit (orb775361)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Porcine |
| Range | 1.57-100 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | P34930 |
| Sensitivity | 0.61 ng/mL |
| Expiration Date | Please enquire. |

