

Product Datasheet

Mouse PKM2 ELISA Kit (orb775380)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Mouse |
| Range | 0.32-20 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma and other biological fluids |
| Uniprot ID | P52480 |
| Sensitivity | 0.124 ng/mL |
| Expiration Date | Please enquire. |

