

Product Datasheet

Gallus IL6 ELISA Kit (orb775471)

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Gallus |
| Range | 15.63-1000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | Q90Y10 |
| Sensitivity | 5.6 pg/mL |
| Expiration Date | Please enquire. |

