

Product Datasheet

Human SPB ELISA Kit (orb775635)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	1.57-100 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, lung lavage fluid, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P07988
Sensitivity	0.61 ng/mL
Expiration Date	Please enquire.

