

## Product Datasheet

### Human NKB ELISA Kit (orb775665)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Human
<b>Range</b>	12.5-800 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<b>Q9UHF0</b>
<b>Sensitivity</b>	3.12 pg/mL
<b>Expiration Date</b>	Please enquire.

