

Product Datasheet

Human NTXI ELISA Kit (orb775707)

| | |
|----------------------------|---|
| Description | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| Reactivity | Human |
| Range | 3.13-200 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Sensitivity | 0.94 ng/mL |
| Expiration Date | Please enquire. |

