

Product Datasheet

Human AGRP ELISA Kit (orb775733)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 7.82-500 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma and other biological fluids. |
| Uniprot ID | 000253 |
| Sensitivity | 3.1 pg/mL |
| Expiration Date | Please enquire. |

