

Product Datasheet

Mouse GHRL ELISA Kit (orb775809)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	156.25-10000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9UBU5
Sensitivity	45.6 pg/mL
Expiration Date	Please enquire.

