

## Product Datasheet

### Gallus TGFb1 ELISA Kit (orb775903)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>          | Gallus  |
| <b>Range</b>               | 15.63-1000 pg/mL  |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Sample Types</b>        | Serum, plasma and other biological fluids.  |
| <b>Uniprot ID</b>          | <b>P09531</b>   |
| <b>Sensitivity</b>         | 6.2 pg/mL   |
| <b>Expiration Date</b>     | Please enquire.   |

