

## Product Datasheet

### Human IL25 ELISA Kit (orb776015)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 15.63-1000 pg/mL  |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Sample Types</b>        | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids    |
| <b>Uniprot ID</b>          | <b>Q9H293</b>   |
| <b>Sensitivity</b>         | 5.9 pg/mL   |
| <b>Expiration Date</b>     | Please enquire.   |

