

Product Datasheet

Mouse bTG ELISA Kit (orb776053)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	1.57-100 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P62325
Sensitivity	0.65 pg/mL
Expiration Date	Please enquire.

