

Product Datasheet

Rat NT3 ELISA Kit (orb776161)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	3.13-200 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P18280
Sensitivity	1.33 pg/mL
Expiration Date	Please enquire.

