

Product Datasheet

Rat EDN1 ELISA Kit (orb776300)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Rat |
| Range | 7.82-500 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | P22387 |
| Sensitivity | 2.71 pg/mL |
| Expiration Date | Please enquire. |

