

## Product Datasheet

### Mouse GHRH ELISA Kit (orb776426)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Mouse
<b>Range</b>	15.63-1000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<b>P16043</b>
<b>Sensitivity</b>	4.63 pg/mL
<b>Expiration Date</b>	Please enquire.

