

## Product Datasheet

### Mouse IL33 ELISA Kit (orb776525)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Mouse
<b>Range</b>	7.82-500 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<b>Q8BVZ5</b>
<b>Sensitivity</b>	3.2 pg/mL
<b>Expiration Date</b>	Please enquire.

