

Product Datasheet

Human APOC2 ELISA Kit (orb776604)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 1.57-100 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma and other biological fluids |
| Uniprot ID | P02655 |
| Sensitivity | 0.74 ng/mL |
| Expiration Date | Please enquire. |

