

Product Datasheet

Gallus IL8 ELISA Kit (orb776621)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Gallus
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P08317
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.

