

Product Datasheet

Gallus IL10 ELISA Kit (orb776644)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Gallus
Range	7.82-500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q6A2H4
Sensitivity	3.3 pg/mL
Expiration Date	Please enquire.

