

## Product Datasheet

### Mouse GC ELISA Kit (orb776651)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Mouse
<b>Range</b>	23.44-1500 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<b>P55095</b>
<b>Sensitivity</b>	7.22 pg/mL
<b>Expiration Date</b>	Please enquire.

