

Product Datasheet

Mouse CFD ELISA Kit (orb776734)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.47-30 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P03953
Sensitivity	0.22 ng/mL
Expiration Date	Please enquire.

