

Product Datasheet

Animal IL8 ELISA Kit (orb776739)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Bovine |
| Range | 15.63-1000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | P79255 |
| Sensitivity | 6.2 pg/mL |
| Expiration Date | Please enquire. |

