

Product Datasheet

Gallus CNTF ELISA Kit (orb776785)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Gallus |
| Range | 62.5-4000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. |
| Uniprot ID | Q02011 |
| Sensitivity | 23.7 pg/mL |
| Expiration Date | Please enquire. |

