

Product Datasheet

Mouse ANGPTL4 ELISA Kit (orb777010)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9Z1P8
Sensitivity	0.055 ng/mL
Expiration Date	Please enquire.

