

Product Datasheet

Human PI3 ELISA Kit (orb777026)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, urine, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P19957
Sensitivity	0.124 ng/mL
Expiration Date	Please enquire.

